

Product Highlights

- MRI Ready device has been tested for safe performance of an MRI scan using a 1,5 T (Tesla) field-strength MRI scanner when used in combination with an MRI Conditional lead^{1,2}
- Improved shape with reduced volume and thickness
- Parylene coating for improved abrasion resistance
- DynamicTx™ Over-Current Detection Algorithm automatically changes shock configurations to ensure delivery of high voltage therapy when high current is detected
- Cold Can programmability provides an additional RV-SVC Shock Configuration to decouple the can from the shocking vector parameters in cases of lead problems
- ShockGuard™ technology with DecisionTx™ programming designed to reduce inappropriate therapy and minimise the need for programming adjustments at implant
 - SecureSense™ RV lead noise discrimination detects sustained and short bursts of lead noise that would otherwise go unnoticed or potentially lead to one or more inappropriate shocks
 - Far Field MD™ morphology discrimination improves SVT and VT discrimination for reduced inappropriate therapies
- Low Frequency Attenuation filter designed to enhance sensing performance and may reduce the possibility of oversensing T-waves
- SenseAbility™ feature provides flexibility to fine-tune programming around T-wave oversensing without decreasing sensitivity
- DF4 connector designed to streamline defibrillation connections into a single terminal pin and reduce the number of set screws
- CorVue™ congestion monitoring feature monitors the intrathoracic impedance in multiple vectors for improved accuracy, and it provides the option for both patient and physician alerts
- Antitachycardia pacing (ATP) while charging and prior to charging in the VF zone further extends the programming options for terminating tachyarrhythmias without a high-voltage shock
- ST monitoring capability provides unprecedented, continuous insight into significant ST shift events and associated ventricular arrhythmias through enhanced monitoring of iEGM and ST-segment as a diagnostic tool to help guide appropriate clinical action
- 36 J delivered energy safety shock option can provide a greater DFT safety margin
- DeFT Response™ technology offers the most noninvasive options for managing high DFTs
- QHR™† chemistry battery provides greater capacity for enhanced longevity and improved charge time performance compared to previous SVO batteries



Merlin@home™
Transmitter
Compatible



Ordering Information

Contents: Single-Chamber Implantable Cardioverter Defibrillator (ICD)

Model Number	Dimensions (H x W x T, mm)	Weight (g)	Volume (cc)	Connector Defibrillation	Connector Sense/Pace
CD1377-36C	68 x 51 x 12	66	31	DF1	IS-1
CD1377-36QC*	66 x 51 x 12	67	30	DF4	DF4

*Indicates models that are MRI Conditional^{1,2}

Indications: The devices are intended to provide ventricular antitachycardia pacing and ventricular defibrillation for automated treatment of life-threatening ventricular arrhythmias.

Contraindications: Contraindications for use of the implantable cardioverter defibrillator include ventricular tachyarrhythmias resulting from transient or correctable factors such as drug toxicity, electrolyte imbalance, or acute myocardial infarction.

Adverse Events: Implantation of the implantable cardioverter defibrillator, like that of any other device, involves risks, some possibly life-threatening. These include but are not limited to the following: acute hemorrhage/bleeding, air emboli, arrhythmia acceleration, cardiac or venous perforation, cardiogenic shock, cyst formation, erosion, exacerbation of heart failure, extrusion, fibrotic tissue growth, fluid accumulation, hematoma formation, histotoxic reactions, infection, keloid formation, myocardial irritability, nerve damage, pneumothorax, thromboemboli, venous occlusion. Other possible adverse effects include mortality due to: component failure,

device-programmer communication failure, lead abrasion, lead dislodgment or poor lead placement, lead fracture, inability to defibrillate, inhibited therapy for a ventricular tachycardia, interruption of function due to electrical or magnetic interference, shunting of energy from defibrillation paddles, system failure due to ionising radiation. Other possible adverse effects include mortality due to inappropriate delivery of therapy caused by: multiple counting of cardiac events including T waves, P waves, or supplemental pacemaker stimuli. Among the psychological effects of device implantation are imagined pulsing, dependency, fear of inappropriate pulsing, and fear of losing pulse capability.

Refer to the User's Manual for detailed indications, contraindications, warnings, precautions and potential adverse events.

† QHR is a trademark of Greatbatch Medical

Ellipse™ VR

Single-Chamber Implantable Cardioverter Defibrillator (ICD)

Product Specifications

Physical Specifications		
Models	CD1377-36C	CD1377-36QC
Telemetry	RF	RF
Delivered/Stored Energy (J)	36/39	36/39
Volume (cc)	31	30
Weight (g)	66	67
Size (mm)	68 x 51 x 12	66 x 51 x 12
Defibrillation Lead Connections	DF1	DF4
Sense/Pace Lead Connections	IS-1	DF4
High-Voltage Can	Electrically active titanium can	Electrically active titanium can
Coating	Parylene	Parylene
MRI Conditional	No	Yes-MRI ready

Parameter	Settings
Sensing/Detection	
SenseAbility™ Technology	Automatic Sensitivity Control adjustment for ventricular events
Low Frequency Attenuation	On; Off
Threshold Start	(Post-Sensed; Ventricular) 50; 62.5; 75; 100%; (Post-Paced; Ventricular) Auto; 0.2-3.0 mV (Post-Sense/Post-Pace; Ventricular) 0-220
Decay Delay	125; 157
Ventricular Sense Refractory (ms)	3 zone programming - 1 zone, 2 zones or 3 zones (VT-1, VT-2, VF)
Detection Zones	Sudden Onset, Interval Stability; Sinus Interval History; Morphology
SVT Discriminators	Discrimination (Far Field MD or Original MD) with Manual (Original MD) or Automatic Template Update
Discrimination modes	On, Passive, Off
SVT Threshold	150-240 min ⁻¹
SVT Timeout	0.25-5 min
Monitor Mode	Detection, discrimination and diagnostics, no therapy delivery (VT or VT-1 zone)
Reconfirmation	Continuous sensing during charging
Lead Noise Discrimination	SecureSense™ RV lead noise discrimination (On; On with Timeout; Passive; Off)

Antitachycardia Pacing Therapy	
ATP Configurations	Ramp; Burst; Scan; 1 or 2 schemes per VT zone
ATP in VF Zone	ATP While Charging; ATP Prior to Charging; Off
ATP Upper Rate Cutoff	150 - 300 min ⁻¹
Burst Cycle Length	Adaptive; Readaptive or Fixed
Min. Burst Cycle Length (ms)	150-400 in increments of 5
Number of Bursts	1-15
Number of Stimuli	2-20
Add Stimuli per Burst	On; Off
ATP Pulse Amplitude (V)	7.5 Independent from Bradycardia and Post-Therapy Pacing
ATP Pulse Width (ms)	1.0 or 1.5 Independently programmable from Bradycardia and Post-Therapy Pacing

High-Voltage Therapy	
DynamicTx™ Algorithm	On; Off
DeFT Response™ Technology	Programmable pulse width for P1/P2 and tilt
High-Voltage Output Mode	Fixed Pulse Width; Fixed Tilt
Waveform	Biphasic; Monophasic
RV Polarity	Cathode (-); Anode (+)
Electrode Configuration	RV to Can; RV to SVC/Can; RV to SVC

Bradycardia Pacing	
Permanent Modes	Off; VVI(R)
Temporary Modes	Off; VVI; VOO
Rate-Adaptive Sensor	On, Off, Passive
Programmable Rate Parameters	Off; Base Rate (min ⁻¹); Rest Rate (min ⁻¹); Maximum Sensor Rate (min ⁻¹); Pulse Amplitude (RV) (V); Pulse Width (RV) (ms); Hysteresis Rate (min ⁻¹); Rate Hysteresis with Search
Ventricular AutoCapture™ Pacing System	On; Off

Post-Therapy Pacing (Independently programmable from Bradycardia and ATP)	
Post-Shock Pacing Mode	Off; VVI
Post-Shock Base Rate (min ⁻¹)	30-100 in increments of 5
Post-Shock Pacing Duration (min)	Off; 0.5; 1; 2.5; 5; 7.5; or 10

Device Testing/Induction Methods	
DC Fibber™ Pulse Duration (sec)	0.5-5.0
Burst Fibber Cycle Length (ms)	20-100
Noninvasive Programmed Stimulation (NIPS)	2-25 stimuli with up to 3 extrastimuli

Customer Support: 46-8-474-4756

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Devices depicted may not be available in all countries. Check with your St. Jude Medical representative for product availability in your country.

Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL and the nine-squares symbol are trademarks and service marks of St. Jude Medical, Inc. and its related companies. © 2015 St. Jude Medical, Inc. All Rights Reserved.

SJM-ELP-1114-0005(1) | Item approved for international use only.

Implantable Cardioverter Defibrillator (ICD) Devices

Patient Notifiers	
Programmable Notifiers (On; Off)	Device at ERI; Charge Time Limit Reached; Possible HV Circuit Damage; Ventricular Lead Impedance Out of Range; High-Voltage Lead Impedance Out of Range; %V pacing; CorVue™ Congestion Trigger; SecureSense lead noise detected, non-sustained lead noise detected, ST Episodes (Type I only)
Device Parameter Reset	On
Entry into Backup VVI Mode	On
Vibration Duration (sec)	2; 4; 6; 8; 10; 12; 14; 16
Number of Vibrations per Notification	2
Number of Notifications	1-16
Time Between Notifications (hours)	10; 22
Electrograms and Diagnostics	
Stored Electrograms	Up to 45 minutes including up to one minute programmable pre-trigger data per VT/VF diagnosis/detection electrograms; triggers include: diagnosis; detection; therapy; PC shock delivery; noise reversion; magnet reversion; and morphology template verification; lead noise detected, non-sustained lead noise detected, NSVT/NSVF
Therapy Summary	Diagram of therapies delivered
Episodes Summary	Directory listing of up to 60 episodes with access to more details including stored electrograms
Lifetime Diagnostics	History of bradycardia events and device-initiated charging
Ventricular HV Lead Impedance Trend	Multi-Vector Trend Data
Histograms	Event Histogram; Ventricular Heart Rate Histogram; Exercise and Activity Trending; DirectTrend™ reports up to 1 year
Real-Time Measurements (RTM)	Pacing lead impedances; high-voltage lead impedances; and signal amplitudes
ST Monitoring	ST Histogram Data; Long-term ST Deviation Trend; ST Episode Log; ST Episode Details; 24-Hour ST & HR Trend; ST EGM Baseline & Snapshots prior to ST Episode, VT/VF, Interrogation (Snapshots and 24-hour trend at time of interrogation)
CorVue™ Congestion Monitoring	On; Off
CorVue Congestion Trigger	8-18 days

MRI Scan Restrictions

Lead Model	Whole Body SAR	Scan Zone Restrictions
Durata™ Lead		If MRI Mode = Pacing OFF: No scan zone restrictions
7120Q (lead lengths: 58 cm, 65 cm)	≤ 2 W/kg	
7122Q (lead length: 58 cm)	≤ 2 W/kg	
7122Q (lead length: 65 cm)	≤ 1,6 W/kg	
Optisure™ Lead		If MRI Mode = VOO or DOO: Isocenter cannot be placed between C1 and L2
LDA220Q (lead lengths: 58 cm, 65 cm)	≤ 2 W/kg	
LDA210Q (lead length: 58 cm)	≤ 2 W/kg	
LDA210Q (lead length: 65 cm)	≤ 1,6 W/kg	

1. MRI Conditional Field Strength: 1,5 Tesla
2. See MRI Procedure Information for approved MR Conditional Systems Device/Lead combinations and scan parameters

